




#### Vorteile unserer Endoskope:

- Kompatibel mit allen Lichtquellen auf dem Markt (universalkompatibel)
- Extrem robuste Rohrkonstruktion mit verstärktem 2-Rohr-Edelstahl-System
- Hochwertig vergütetes, hochauflösendes Linsen-System
- Optische Saphirglas-Frontlinse, kratzfest und hitzebeständig fixiert mit Epoxidharz
- Hitzebeständiges Polymer-Okular (PEEK)
- Hohe Dichtigkeit durch gelaserte Schnittstellen
- Mehr Bildinformation durch großen Öffnungswinkel
- Hohe Bildschärfe im gesamten Bildbereich
- Gleichmäßige Helligkeit bis zum Rand durch Lichtkonus
- Autoklavierbar
- 24 Monate Garantie

#### Advantages of our endoscopes:

- Compatible to all light sources available (universal camera compatible)
- Extremely robust stainless steel tube construction with reinforced outer tube
- High precision quality lens system
- Optical sapphire glass front lens, scratch resistant and fixed with heat-resistant epoxide
- Heat-resistant polymer ocular (PEEK)
- Tight sealing due to laser welded joints
- More image information due to wide viewing angle
- High focus quality in entire viewing area
- Constant brightness up to the edge of viewing area
- Autoclavable
- 24 months guarantee



**All unsere Endoskope sind aus Phynoxstahl**  
für verbesserten Schutz gegen Verbiegen und Bruch

**All our endoscopes are made of phynox steel**  
for improved protection against bending or breaking


## 02 Ohr-Endoskope

### Ear Endoscopes




**Otoskop, HD, autoklavierbar / Otoscope, HD, autoclavable**

| Art.-Nr.<br>Art.-No. | Größe<br>Size (mm) | Länge<br>Length (cm) | Winkel<br>Angle |
|----------------------|--------------------|----------------------|-----------------|
| 10-102-00            | 2,7                | 5,0                  | 0°              |
| 10-102-30            | 2,7                | 5,0                  | 30°             |
| 10-104-00            | 2,4                | 9,8                  | 0°              |
| 10-104-30            | 2,4                | 9,8                  | 30°             |




**Gehörgangs-Endoskop, HD, autoklavierbar, AD 4 mm, 5 cm  
Ear canal endoscope, HD, autoclavable, OD 4 mm, 5 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 10-106-00    | 0°             |
| 10-107-30    | 30°            |




**Gehörgangs-Endoskop mit Ohrspekulum,  
HD, autoklavierbar, AD 4 mm, 5 cm  
Ear canal endoscope with ear speculum,  
HD, autoclavable, OD 4 mm, 5 cm**

| Art.-Nr./No. | Winkel / Angle                                               |
|--------------|--------------------------------------------------------------|
| 10-108-00    | 0°                                                           |
| 10-108-99    | Ersatzspekulum für 10-108-00<br>spare speculum for 10-108-00 |



**Otoskop, HD, Großbild, autoklavierbar, AD 3 mm, 11 cm  
Otoscope, HD, large picture, autoclavable, OD 3 mm, 11 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 10-109-00    | 0°             |
| 10-109-30    | 30°            |
| 10-109-45    | 45°            |
| 10-109-70    | 70°            |




**Otoskop, HD, Großbild, autoklavierbar, AD 3 mm, 14,1 cm  
Otoscope, HD, large picture, autoclavable, OD 3 mm, 14.1 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 10-090-00    | 0°             |
| 10-090-30    | 30°            |
| 10-090-45    | 45°            |

## 03 NNH-Endoskope


### Sinus Endoscopes



**Otoskop, HD, Großbild, autoklavierbar, AD 2,4 mm, 11 cm**

Otoscope, HD, large picture, autoclavable, OD 2.4 mm, 11 cm


| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 10-111-00    | 0°             |
| 10-111-30    | 30°            |
| 10-111-45    | 45°            |
| 10-111-70    | 70°            |



**Sinuskop, HD, semiflex, autoklavierbar, AD 2,7 mm, 11 cm**

Sinuscope, HD, semiflex, autoclavable, OD 2.7 mm, 11 cm

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-135-00    | 0°             |
| 30-135-30    | 30°            |



**Sinuskop, HD, starr, autoklavierbar, AD 4,0 mm, 17,5 cm**

Sinuscope, HD, rigid, autoclavable, OD 4.0 mm, 17.5 cm

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-142-45    | 45°            |
| 30-142-90    | 90°            |



**Sinuskop, HD 4K, autoklavierbar, AD 4,0 mm, 17,5 cm**


Sinuscope, HD 4K, autoclavable, OD 4.0 mm, 17.5 cm

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-143-00    | 0°             |
| 30-143-30    | 30°            |
| 30-143-70    | 70°            |



### 03 NNH-Endoskope


#### Sinus Endoscopes



**Sinuskop, HD, starr, autoklavierbar, AD 2,7 mm, 18,75 cm**

**Sinuscope, HD, rigid, autoclavable, OD 2.7 mm, 18.75 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-130-00    | 0°             |
| 30-130-30    | 30°            |
| 30-130-70    | 70°            |



**Sinuskop, HD, starr, autoklavierbar, AD 2,7 mm, 14,1 cm**

**Sinuscope, HD, rigid, autoclavable, OD 2.7 mm, 14.1 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-132-00    | 0°             |
| 30-132-30    | 30°            |
| 30-132-45    | 45°            |



**Sinuskop universal, HD, starr, autoklavierbar, AD 2,7 mm, 11 cm**

**Universal sinuscope, HD, rigid, autoclavable, OD 2.7 mm, 11 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 30-134-00    | 0°             |
| 30-134-30    | 30°            |
| 30-134-45    | 45°            |
| 30-134-70    | 70°            |

## 04 Laryngoskopie-Endoskope

### Laryngoscopy Endoscopes



**Laryngoskop, HD, autoklavierbar, AD 4,0 mm, 17,5 cm**

**Laryngoscope, HD, autoclavable, OD 4.0 mm, 17.5 cm**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 40-100-70    | 70°            |
| 40-100-90    | 90°            |



**Laryngoskop, HD, autoklavierbar, AD 10 mm, 19 cm, Handgriff**


**Laryngoscope, HD, autoclavable, OD 10 mm, 19 cm, handle**

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 40-103-70    | 70°            |
| 40-103-90    | 90°            |



## 04 Laryngoskopie-Endoskope

### Laryngoscopy Endoscopes



**Zoom Laryngoskop, HD, einlegbar, AD 10 mm, 20 cm, Handgriff**

Zoom Laryngoscope, HD, soakable, OD 10 mm, 20 cm, handle


| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 40-105-70    | 70°            |
| 40-105-90    | 90°            |



**Laryngoskop, HD, autoklavierbar, AD 8,0 mm, 19,8 cm, 70°**

Laryngoscope, HD, autoclavable, OD 8.0 mm, 19.8 cm, 70°


40-107-70



**Laryngoskop, HD, autoklavierbar, AD 2,7 mm, 30,2 cm**

Laryngoscope, HD, autoclavable, OD 2.7 mm, 30.2 cm

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 40-110-00    | 0°             |
| 40-110-30    | 30°            |
| 40-110-70    | 70°            |



**Laryngoskop, HD, autoklavierbar, AD 4,0 mm, 30,2 cm**

Laryngoscope, HD, autoclavable, OD 4.0 mm, 30.2 cm

| Art.-Nr./No. | Winkel / Angle |
|--------------|----------------|
| 40-112-00    | 0°             |
| 40-112-30    | 30°            |
| 40-112-70    | 70°            |